# Rotor F-55-16-5-PCR – Gebrauchsanweisung

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### 1 Introduction

Rotor F-55-16-5-PCR is a fixed-angle rotor for centrifuging 16 PCR tubes (0.2 ml) or two 5x or 8x rows of 0.2 ml PCR tubes in the Mini Spin<sup>®</sup>, MiniSpin<sup>®</sup> plus and MiniSpin<sup>®</sup> SPACE. The field of application for this rotor is the research, routine and training laboratory in the fields of biosciences, medicine, clinical diagnostics, industry and chemistry.

**Before using the rotor for the first time, please read this information sheet.**

### 2 Safety information

For your personal safety, please be sure to comply with the following regulations unconditionally.

- Rotor F-55-16-5-PCR must only be used in the MiniSpin<sup>®</sup>, MiniSpin<sup>®</sup> plus and MiniSpin<sup>®</sup> SPACE centrifuges for the fields of application mentioned. Use only original Eppendorf spare parts, rotors and rotor lids.
- Rotors and rotor lids must always be properly secured. The centrifuge may only be operated with the rotor firmly tightened.
- Protect the rotor and the rotor lid from mechanical damage. Even minor scratches and cracks can result in serious internal material damage. Avoid damage from aggressive chemicals including, among other things, strong and weak alkalis, strong acids, solutions of mercury, copper and other heavy metal ions, chlorinated hydrocarbons, concentrated salt solutions and phenol. If the rotor is contaminated by aggressive substances, clean it immediately with a neutral agent (e.g. Extran<sup>®</sup> neutral, RBS<sup>®</sup> neutral or Teepol<sup>®</sup> 610 S).
- Do not use rotors and rotor lids with clear signs of corrosion or mechanical damage. You should accordingly check the rotor and the rotor lid regularly.
- Seal tube lids down properly before centrifuging. Open tube lids may rip off during centrifugation and damage the centrifuge, the rotor and the rotor lid.
- When handling toxic liquids or pathogenic microorganisms of risk group II (see World Health Organization: Laboratory Biosafety Manual), comply with the relevant national regulations. If the named liquids are spilled in the rotor or rotor chamber, the centrifuge must be thoroughly and professionally cleaned. Before any cleaning or disinfecting method other than that in Item 5 "Cleaning" is used, please check with Eppendorf that the intended method will not damage the rotor and the rotor lid.

### 3 Delivery package

1 Rotor F-55-16-5-PCR incl. rotor lid and rotor nut (see also Item 4 "Operation")

# Rotor F-55-16-5-PCR – Instructions for use

## 4 Operation

### Fitting and removing the rotor

Place the rotor on the motor shaft followed by the rotor nut. Tighten the rotor nut by turning it clockwise. Before each start, check that the rotor is firmly tightened. To undo the rotor: turn the rotor nut anticlockwise.

### Loading the rotor

The rotor must always be loaded symmetrically. Minimize differences in weight between the filled sample tubes - taring with a scale is recommended. This will reduce wear on the drive and cut running noise.

16 PCR tubes (0.2 ml) or two 5x or 8x row 0.2 ml PCR tubes with the following performance data can be centrifuged in the rotor:

	Relative centrifugal force (= rcf)	Speed in RPM	Max. centrifugal radius (cm)
MiniSpin®	9,840	13,400	4.9
MiniSpin® plus / SPACE	11,520	14,500	4.9

\* rcf = Relative centrifugal force

Maximum load (tubes and contents) per 8x row with this rotor is **3.5 g**.

The rotor may only be centrifuged with the rotor lid tightened in accordance with specification. To do this, push the rotor lid with the seal onto the rotor nut so that the edge of the rotor lid is located in the circumferential groove of the rotor.

## 5 Cleaning

The rotor and the rotor lid can be cleaned with a neutral agent (e.g. Extran® neutral, RBS® neutral, Teepol® 610 S). Bacillo® AF, Meliseptol® und Perform® are recommended for cleaning and disinfecting the rotor and the rotor bores.

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Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA