## Thermal Printer – Instructions for use

### 1 Product description

The thermal printer has one parallel and one serial interface (Centronics and RS 232 C).

Fig. 1: Thermal Printer

### 2 Using this manual

- ▶ Please read this manual completely before using the device for the first time.
- ▶ Please view this manual as part of the product and keep it somewhere easily accessible.
- ▶ When passing the device on to third parties, be sure to include this operating manual.
- ▶ If this manual is lost, please request another one. The current version can be found on our website [www.eppendorf.com](http://www.eppendorf.com) (International).

### 3 Setting DIP switches

- ❗ The DIP switches are preset for use with the following devices: BioPhotometer, BioPhotometer plus, BioSpectrometer basic, BioSpectrometer kinetic. For details, refer to the BioPhotometer/BioSpectrometer's instructions for use.

**New settings** can be entered into the printer program.

- ▶ **POWER** off – press and hold down the **ON LINE** key – **POWER** on.
- ▶ The printout of the current settings starts – release the **ON LINE** key.
- ▶ After the printout has finished, press the **ON LINE** key to start the setting procedure.
- ▶ All settings (new and old) must be re-entered.
- ▶ The settings are switched on by pressing the **ON LINE** key and switched off by pressing the **FEED** key.
- ▶ At the end of entries for each of the three setting ranges, the entry is continued by pressing the **ON LINE** key or completed by pressing the **FEED** key.

**Then switch printer and device off and on again.**

## Thermal Printer – Instructions for use

### 4 Disposal

Always observe the statutory requirements when disposing of the product.

#### Information on the disposal of electrical and electronic devices in the European Community:

Within the European Community, the disposal of electrical devices is regulated by national regulations based on EU Directive 2002/96/EC pertaining to waste electrical and electronic equipment (WEEE).

According to these regulations, any devices supplied after August 13, 2005, to the business-to-business area this product is assigned to, may no longer be disposed of in municipal or household waste. To document this, they have been marked with the following identification:

Because disposal regulations may differ from one country to another within the EU, please contact your supplier if necessary.

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Eppendorf AG · Hamburg · Germany · E-Mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)