Additional Infos (Thermomixer)

11.01.2008

Improved spare part for Thermomixer comfort 5355:

With immediate effect a new silicone ring to fix the heating / cooling plate connections is installed. The order number is 5355 866.060. The formerly used O-ring with part no. 0013 020.132 is no longer available. Please update your spare parts list accordingly.

15.03.2002

There is an operational mode with the new switching power supplies (SIM P069), which may cause some irritation. If the mixer is switched off and immediately switched on again, the beeper may sound several times. After a short time, the device will operate normally.

07.02.00

Thermomixer 5350 compact is emitting a beeping noise?

How to solve.

- Switch on device.
- Set MIX = 500 rpm.
- Switch off device.
- Switch ON device while pressing **Temp-down** key until message SERVICE ON! is displayed.
- Press **Temp-up** key until 6. Heat PWM is displayed.
- Press **Start/Stop** key.
- Press **Temp-up** key to select parameter Hc: PWM Mode.
- Press **Mix-down** key for 5 seconds to switch from FAST mode to SLOW mode.
- Make sure that SLOW mode is active.
- Switch off device.
- Switch on and check device.

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