


Contact person at HQ, Hamburg:	Marco Carmon
Phone:	+49 40 53801-572
E-mail:	Carmon.m@eppendorf.de

## MixMate

### MixMate as of S/N 15001



Vortex mats starting with S/N 10001 are glued to the plate holder to prevent liquid from entering the device. That meant that the vortex mat could no longer be removed and the screws below it, which are used to disassemble the device, could not be removed. Repair or service of the MixMate was temporarily unavailable.

A modified plate holder will be used starting with S/N 15001. The glued vortex mat now has cover caps, which can be removed as necessary, on the screw positions (see figure), which allow the MixMate to be disassembled or repaired. Devices which were manufactured before this change (S/N 10001 to 15000) can now be repaired again. The vortex mat must be cut open on the screw positions to disassemble the plate holder. Therefore, the old plate holder cannot continue being used; it must be replaced by the new plate holder with glued vortex mat and cover caps.

The PCB MixMate has also been changed. The new PCB MixMate is backwards compatible and can be used in all MixMates starting with S/N 10001 without any limitations.

Starting with S/N 10001, a new screening plate is used above the LP MixMate, to replace the old cover of the PCB.

These changes have also brought about changes to the MixMate spare parts list:

Omitted spare parts/invalid order numbers:

Vortex mat	5353 863.004
Plate holder	5353 862.008
Cover PCB MixMate	5353 850.140
PCB MixMate 120 V	5353 855.109
PCB MixMate 240 V	5353 854.102

New spare parts/new order numbers:

Plate holder (incl. vortex mat and cover caps)	5353 862.105
Cover caps (sets of 4)	5353 863.101
Screening plate	5353 850.107
PCB MixMate 120 V	5353 855.206
PCB MixMate 240 V	5353 854.200