

Service Check List

Mastercycler® Essential Check



epServices
for premium performance

Instrument:	Instrument nick name:	Serial number:
Firmware version: Additional instrument software:	Installation date:	Location:
Job number:	Contract number:	Contact person:

VISUAL INSPECTION

Checking exterior	passed	failed
Checking accessories	passed	failed

FUNCTIONAL CHECK (with service functions)

Display	n.a.	passed	failed
Status LED	n.a.	passed	failed
Keyboard	n.a.	passed	failed
Thermoblock	n.a.	passed	failed
Heated Lid	n.a.	passed	failed
<i>realplex</i> module	n.a.	passed	failed
Fan (s)	n.a.	passed	failed
Final test run	---	passed	failed

COMMENTS

Completed service sticker attached to instrument case	YES	NO	

Date completed:	Service Engineer:	Service Engineer Signature:
	Customer:	Customer Signature:

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · 22331 Hamburg · Germany
 E-mail: eppendorf@eppendorf.com

Eppendorf®, the Eppendorf logo, Mastercycler® and epServices® logo are registered trademarks of Eppendorf AG, Germany.
 U.S. Design Patents are listed on www.eppendorf.com/ip. All rights reserved incl. graphics and pictures. GB1/1215. Copyright 2015 © by Eppendorf AG.