



Centrifuge 5417 R – Installation

These instructions supplement the Centrifuge 5417 R operating manual. They do not replace it. Be sure to read the operating manual before you begin to use the centrifuge. The current operating manual version is available on the Internet at www.eppendorf.com.


1.




4.




2.




5.



3.



6.



eppendorf

Centrifuge 5417 R – Installation

Remove transport safety device!

* After unpacking, wait 4 hours before switching on the device.

Transportsicherung entfernen!

* Gerät nach dem Auspacken 4 Stunden still stehen lassen.

Retirer le support d'expédition!

* Après déballage, ne pas déplacer ou mettre en marche l'unité pendant au moins 4 heures.

Rimuovere il dispositivo di fissaggio per il trasporto!

* L' attrezzatura deve restare ferma per 4 ore dopo essere stata disimballata.

¡Desmante el soporte de transporte!

* Después de desembalar, no mueva ni conecte la unidad durante, por lo menos, 4 horas.

Verwijder de transportriem!

* Het apparaat, na het uitpakken, gedurende ten minste 4 uur niet verplaatsen of aanzetten.

Fjern transportafstivningen!

* Enheden må ikke flyttes eller startes før tidligst 4 timer efter udpakning.

Poista pakkausvanteet!

* Älä liikuta tai yritä käynnistää laitetta ennen kuin neljä tuntia on kuunut pakkauksen avaamisesta.

Απομακρύνετε την συσκευασία.

* Μετά την απομάκρυνση της συσκευασίας μην μετακινήσετε και μην θέσετε σε λειτουργία την συσκευή για τουλάχιστον 4 ώρες.

Retire o sistema de imobilização para transporte.

* Depois de desempacotar deixe a unidade imóvel durante pelo menos 4 horas antes de a ligar.

Avlägsna transportstödet!

* Flytta inte enheten eller slå på den på åtminstone fyra timmar efter uppackningen.

eppendorf

Your local distributor: www.eppendorf.com/worldwide

Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: eppendorf@eppendorf.com

Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA

Tel: +1 516 334 7500 · Toll free phone: +1 800-645-3050 · Fax: +1 516 334 7506 · E-mail: info@eppendorf.com

Application Support Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)
support@eppendorf.com

North America: Tel: +1 800 645 3050 · E-mail: techserv@eppendorf.com

Asia Pacific: Tel: +60 3 8023 6869 · E-mail: support_asiapacific@eppendorf.com