


Contact person in HQ, Hamburg:	Ralf Teichmann
Phone:	+494053801-667
E-mail:	teichmann.r@eppendorf.de

## New software version 3.5.0.0 for the BioSpectrometer family 613x and BioPhotometer D30 6133 available

A new software version for the BioSpectrometer family 613x and the BioPhotometer D30 6133 is available for download from [www.eppendorf-support.com](http://www.eppendorf-support.com).


On some devices the self-test did not work correctly, no directory was available in the "Method Selection" screen and no predefined methods were visible. See the following image.



These errors have been corrected in the new software version 3.5.0.0. The installation instructions are also available at [www.eppendorf-support.com](http://www.eppendorf-support.com).

## Neue Software Version 3.5.0.0 für die BioSpectrometer Familie 613x und BioPhotometer D30 6133 verfügbar

Für die BioSpectrometer Familie 613x und das BioPhotometer D30 6133, steht eine neue Software Version zum Herunterladen auf der [www.eppendorf-support.com](http://www.eppendorf-support.com) zur Verfügung. Bei einigen Geräten kam es zu einem fehlerhaften Selbsttest und im Bildschirm „Method Selection“, war kein Verzeichnis vorhanden und die vordefinierten Methoden waren ebenfalls nicht sichtbar. Siehe Bild.



Mit dieser Software Version 3.5.0.0 wird dieser Fehler beseitigt. Die Anweisung zum Installieren der Software finden Sie ebenfalls auf der [www.eppendorf-support.com](http://www.eppendorf-support.com).