

| | |
|--------------------------------|--------------------------|
| Contact person in HQ, Hamburg: | Ralf Teichmann |
| Phone: | +494053801-667 |
| Email: | teichmann.r@eppendorf.de |

A new Ethernet-enabled mainboard with a new software version 4.2.3.0 for the BioSpectrometer family 613x

A new “PCB BioSpectrometer” with an Ethernet interface is available for the devices in the BioSpectrometer family 613x. The order no. for this Ethernet-enabled PCB is: 6136 859.300.

This new PCB is already used in devices as of the following serial number:

BioSpectrometer „basic” Seriennummer: 6135xx702304
BioSpectrometer „kinetic” Seriennummer: 6136xx702304
BioSpectrometer „fluorescence” Seriennummer: 6137xx702304

No Ethernet interface is available for the BioPhotometer D30. Use the PCB BioSpectrometer with the order no. 6136 859.003 to repair the BioPhotometer D30 and the older BioSpectrometers up to the serial number

BioSpectrometer „basic” Seriennummer: 6135xx702303
BioSpectrometer „kinetic” Seriennummer: 6136xx702303
BioSpectrometer „fluorescence” Seriennummer: 6137xx702303

A backward-compatible software (version 4.2.3.0.) will be made available for both versions of the “PCB BioSpectrometer”.

The software can be downloaded from eppendorf-support.com.