

Contact person in HQ, Hamburg:	Ingo Dietrich
Phone:	+ 49 40 - 538 01 - 314
Email:	dietrich.i@eppendorf.de

## New software version for Control Panel 6320

Starting with lot no. A242767G (serial number 6320AG004941), all Control Panels 6320 carry the new software version 1.014.

The software version for the Mastercycler pro and Mastercycler ep thermo module remains the same (version 3.024).

The changes in version 1.014 are:

- adaption to a 32 bit flash memory drive
- change of the acoustic signal frequency (1.2 Hz to 2.4 Hz)

The software version is backward compatible with all Control Panels 6320.

---

## Neue Software-Version für Control Panel 6320

Ab Lotnummer A242767G (Seriennummer 6320AG004941) werden alle Control Panel 6320 mit der neuen Software-Version 1.014 ausgeliefert.

Die Software für die Mastercycler pro- und Mastercycler ep- Thermomodule ändert sich nicht (Version 3.024 bleibt bestehen).

Die Änderungen in Software-Version 1.014 sind:

- notwendige Anpassung an einen neuen 32 MB Flash-Speicher
- Änderung der Pieperfrequenz (Verdopplung auf von 1,2 Hz auf 2,4 kHz)

Die neue Software-Version ist rückwärtskompatibel zu allen Control Panels 6320.