

Contact person in HQ, Hamburg:	Uwe Krohss
Phone:	+49 40-53801-193
Email:	Krohss.u@eppendorf.de

For the expanded Mastercycler nexus family, you can find on www.eppendorf-support.com:

- software V 2.2.0.0. including error listing V 2.2.0.0. for download
- updated service manual V 2.2.0.0 for following Mastercycler:

Model	Mastercycler	Thermoblock
6330	MC nexus flat eco	Flat, without control panel
6331	MC nexus gradient	96 Well gradient Alu
6332	MC nexus eco	96 Well Alu, without control panel
6333	MC nexus	96 Well Alu
6334	MC nexus eco gradient	96 Well gradient Alu, without control panel
6335	MC nexus flat	Flat
6345	Mastercycler nexus GSX1	96 Well gradient Silver
6346	Mastercycler nexus SX1	96 Well Silver
6347	Mastercycler nexus GSX1e	96 Well gradient Silver, without control panel
6348	Mastercycler nexus SX1e	96 Well Silver, without control panel

Verification of the Mastercycler nexus flat is only possible with the Eppendorf Multi-Channel-Verification System and sensor for Flat Thermoblocks.