

Contact person in HQ, Hamburg:	Uwe Krohss
Phone:	+494053801-193
Email:	Krohss.u@eppendorf.de

Improved Thermoblocks for Mastercycler nexus are available

In the case of errors 412 - 414 and 216 – 220 we recommend to insert new Thermoblocks with following order numbers.

Order no. for

96 AI Thermoblock 6331 864.003 was replaced by **6331 864.500**,

96 Ag Thermoblock 6345 864.006 was replaced by **6345 864.502**,

Flat – Thermoblock 6331 863.007 was replaced by **6331 863.503** and

X2 – Thermoblock 6336 864.001 was replaced by **6336 864.508**.

Please change your order information and download the updated Service Manual from Support.com.

Verbesserte Thermoblöcke für die Mastercycler nexus sind verfügbar

Wir empfehlen wenn die Fehler 412 - 414 und 216 – 220 auftreten, neue Thermoblöcke mit folgenden

Bestellnummern einzusetzen.

Die Bestellnummer für den,

96 AI Thermoblock 6331 864.003 wurde ersetzt in **6331 864.500**,

96 Ag Thermoblock 6345 864.006 wurde ersetzt in **6345 864.502**,

Flach – Thermoblock 6331 863.007 wurde ersetzt in **6331 863.503** und

X2 – Thermoblock 6336 864.001 wurde ersetzt in **6336 864.508**.

Bitte ändern Sie Ihre Bestellinformationen und laden sich von der Support.com das geänderte Servicemanual.