

Contact person in HQ, Hamburg:	Uwe Krohss
Phone:	+494053801-193
Email:	Krohss.u@eppendorf.de

Mastercycler nexus service manual

V 3.5.1.0. is available on www.eppendorf-support.com

Updates are:

Disassembly for PCB Display, PCB Control pad, Thermoblock revised.

Error search for errors 1441-1444 revised.

Spare parts Heating frame, E-PAC parts, Switched-mode power supply revised.

New software version.

Note: A new Switched-mode power supply will be installed in Mastercycler nexus in 2016. This new Switched-mode power supply 6331 898.005 will be available as a service spare part In 2017. The Switched-mode power supply 6331 860.008 is still available.

More information in the new Service manual in section 6.7.29 Mounting the PCB Switched-mode mains/power supply.

Mastercycler nexus Servicemanual

V 3.5.1.0. ist verfügbar auf www.eppendorf-support.com

Was überarbeitet wurde:

Demontage für LP Anzeige, LP Bedienteil, Thermoblock überarbeitet.

Fehlersuche für Fehler 1441-1444 überarbeitet.

Ersatzteile Heizrahmen, E-PAC-Teile, Schaltnetzteil überarbeitet.

Neue Software-Version.

Hinweis: Ein neues Schaltnetzteil wird in den Mastercycler nexus in 2016 verbaut.

Diese neue Schaltnetzteil 6331 898.005 wird als Service Ersatzteil erst in 2017 verfügbar sein. Das Schaltnetzteil 6331 860.008 ist weiter verfügbar.

Weitere Informationen im neuen Servicemanual im Kapitel 6.7.29 LP Schaltnetzteil montieren.