

| | |
|--------------------------------|-------------------------|
| Contact person in HQ, Hamburg: | Ingo Dietrich |
| Telefon: | + 49 40 - 538 01 - 314 |
| Email: | dietrich.i@eppendorf.de |

A new PCB Power supply unit (6321 856.509) is available for Mastercycler pro 6321 / 6325 Mastercycler ep 5341 / 5344 / 5345 and Mastercycler realplex 6300 / 6302.

The new PCB Power supply unit (6321 856.509) replaces the previous PCB Power supply 6321 856.002 and 5341 856.002.

We recommend that you include this information together with Chapter 2 of your Mastercycler Service Manual, to have all relevant product information available in one place.

Für den Mastercycler pro 6321 / 6325, Mastercycler ep 5341 / 5344 / 5345 und den Mastercycler realplex 6300 / 6302 ist eine neue LP-Netzteil (6321 856.509) verfügbar.

Die neue LP-Netzteil (6321 856.509) ersetzt die bisherigen LP-Netzteil 6321 856.002 und 5341 856.002 .

Wir empfehlen, diese Info Ihrem Service Manual Mastercycler Kapitel 2 beizufügen, damit alle relevanten Produktinformationen zusammenhängend abgelegt sind.