

We recommend that you include this information together with Chapter 12 of your Service Manual for the Thermomixer/Thermostat, so that all relevant product information can be found in one place.

## **Thermomixer compact 5350 Thermomixer comfort 5355**

### **Oscillating frame system compact & comfort**

The spare parts oscillating frame system compact and oscillating frame system comfort (see order numbers below) are no longer manufactured. Both spare parts can no longer be ordered after the sale of existing inventories. A new spare part, which can be used for both devices, has been developed with a new order number in order to serve as a replacement for both oscillating frame systems. The new spare part "Oscillating frame 5350/5355" can be used in all Thermomixer compact and Thermomixer comfort units without limitation.

Order numbers **old** (after sale of existing inventories no longer available):

5350 867.009      Oscillating frame system compact  
5355 866.000      Oscillating frame system comfort

Order numbers **new**:

5350 867.106      Oscillating frame 5350/5355

---

Wir empfehlen, diese Info Ihrem Service Manual Thermomixer/Thermostat, Kapitel 12, beizufügen, damit alle relevanten Produktinformationen zusammenhängend abgelegt sind.

## **Thermomixer compact 5350 Thermomixer comfort 5355**

### **Schwingrahmensystem compact & comfort**

Die Ersatzteile Schwingrahmensystem compact und Schwingrahmensystem comfort (Bestellnummern siehe unten) werden nicht mehr produziert. Nach Abverkauf der Lagerbestände werden beide Ersatzteile nicht mehr bestellbar sein.

Als Ersatz für beide Schwingrahmensysteme wurde ein für beide Geräte einsetzbares Ersatzteil mit neuer Bestellnummer geschaffen. Das neue Ersatzteil „Schwingrahmen 5350/5355“ kann ohne Einschränkung in allen Thermomixer compact und comfort eingesetzt werden.

Bestellnummern **alt** (nach Abverkauf nicht mehr bestellbar):

5350 867.009      Schwingrahmensystem compact  
5355 866.000      Schwingrahmensystem comfort

Bestellnummer **neu**:

5350 867.106      Schwingrahmen 5350/5355